



Molecular Imaging Products Company

Division of Summit Anesthesia Support

ANESTHESIA TECHNOLOGIES

ProSTATION[®]

FOR SMALL ANIMAL SURGERY



AS-01-0491

ProStation Kit
AS-01-0491-K

[CLICK HERE FOR MORE INFORMATION](#)

- ❖ **Compact 5" wide x 13" length x 3/4" high**
- ❖ **Heated work surface**
- ❖ **Downdraft anesthetic gas evacuation**
- ❖ **Autoclaveable, non-reactive to most disinfectants**
- ❖ **Works with all non-rebreathing systems**

63043 Lower Meadow Drive, Unit #110 | Bend, OR 97701 | toll free 800-877-8989 | phone 541-385-9967 | fax 541-385-9662

www.mipcompany.com